## VACUMATIC® Vacuum Technology

With Vacumatic vacuum technology, you can guarantee hygienic waste disposal, while facilitating a pleasant living and working environment.



### **The Vacumatic Technology**

Vacumatic technology removes the air from waste bags and then hermetically seals them removing the problem of exposed bacteria, viruses and unpleasant odors. The Vacumatic system provides a hygienic, easy to operate solution to bad odor in health care facilities.

The patented bags that are used in the Vacumatic system have a strong wall structure which prevents leakage of air or liquid after sealing. They are available in clear or yellow & black 'tiger' bags.

The Vacumatic system aids collection, storage and transportation of incontinence and other clinical waste in a safe, easy to operate and economical way.

If your facility relies on incontinence products, Vacumatic allows you to dispose of waste quickly, easily and hygienically.

The Vacumatic system comprises of multiple collection trolleys and a central vacuum unit which once each special bag is full, it will be compressed and heat sealed to prevent spillage or odor escape.



#### **Benefits**

- Hands free operation reducing the risk of cross-contamination
- Suitable for use in combination with all types of incontinence products
- Unit is fitted with bacteria and odor filters for the protection of staff and the environment
- The strong Vacumatic bags negate the risk of splitting during movement to outside garbage areas and reduces odor whilst waiting for collection
- Plugs into standard electric wall socket







# VACUMATIC® Vacuum Technology

## **VACUMATIC®**

**Vacuum Technology** 



## **Vacumatic Trolley**

The Vacumatic Trolley has been specifically designed as part of the Vacumatic waste disposal system.

The Vacumatic Trolley uses a unique rubber clamping system operated by a foot pedal which opens the clamp spring top and closes automatically when the foot is removed. This means that the waste bag is never left open, preventing any bad odors from escaping, keeping your facility completely safe, efficient and odor free.

DDC Dolphin Order Code	Description	Qty
VAC0007	Vacumatic Trolley	1





## **Vacumatic Bags**

Vacumatic bags are strong (0.004" thick material) and have a patented structure to maintain strength. They will not allow the escape of liquids or foul odor once sealed.

DDC Dolphin Order Code	Description	Qty
VAC0009	Vacumatic Trolley Bags - Clear	100
VAC0010	Vacumatic Trolley Bags - Tiger	70









VAC0001

Vacumatic

## VACUMATIC® Vacuum Technology

#### **Vacumatic Starter Pack**

With the Vacumatic system, you can guarantee hygienic waste disposal, whilst facilitating a pleasant living and working environment with no more bad smell from exposed waste.

The Vacumatic Starter Pack contains all that is needed to start improving hygiene standards and reducing bad odor.

The pack consists of

- 1 x Vacumatic waste disposal unit.
- 2 x Vacumatic collection trolleys.
- 4 cartons of Yellow & Black tiger bags (3 months supply based on three bags per day).

DDC Dolphin Order Code	Description	Qty
VACUSA3M	1 x Vacumatic waste disposal unit, 2 x Vacumatic collection trolleys., 4 cartons of Yellow & Black tiger bags.	1

### **Machine Consumables**

The Vacumatic machine requires an annual maintenance visit to replace wear items and ensure that it continues to offer effective hygiene management.



DDC Dolphin Order Code	Description	Qty
VACKIT12	12 month service pack - contains replacement odor filter and foam sealing strip	1
VACKIT24	24 month service pack - contains replacement odor filter, bacteria filter, foam sealing strip and heat sealing bar	1

This document has been carefully formulated to be a true representation of the product, however the manufacturer accepts no responsibility for any changes, or for any consequence arising from this document.



