



# Bedpan Washer Consumables



## EcoWash+ Wash and Scale Inhibitor

**By far the biggest enemy of bedpan washer disinfectors is limescale, especially in hard water areas.**

Limescale can build up in the steam generator (notably on the heating element) and in the pipes, reducing the flow of water. This build-up of limescale will cause a bedpan washer to work less efficiently and require longer cycles – wasting energy – until it eventually breaks down. To continue to run efficiently, a bedpan washer needs limescale inhibitor.

### Hygenex EcoWash+ 5L



#### Product Description

Wash and scale inhibitor specially formulated for use in DDC Dolphin Panamatic washer disinfectors. EcoWash+ supports the daily functioning of washer disinfectors and helps prolong the life of machinery through the support of scale management.

Order Code	Description	Pack Qty
99-01-052	EcoWash+ Wash & Scale Inhibitor 5L	2

**The Scalematic limescale inhibitor reduces the build-up of limescale throughout your bedpan washer.**

The Scalematic cartridge is fitted on the water supply prior to the machine connection point. By introducing polyphosphate into the water feed and sequestering the calcium in the water it helps prevent the formation of scale.

Replacing your Scalematic every 6 months will ensure your machinery stays in good working order for longer.

### Scalematic



#### Product Description

The Scalematic limescale inhibitor safeguards your machinery against limescale build-up.

Order Code	Description	Qty
SCALEMATIC-SYS	Scalematic Limescale Inhibitor Unit	1
90-01-003	Scalematic Limescale Inhibitor replacement cartridge	1

## Chemical Cabinets

**Our wall mounted chemical cabinets are made of powder coated steel and are available with standard key lock or combination lock.**

The cabinet is supplied with adaptor and flexible conduit to support the chemical supply tube linked directly to the machine.

### Lockable Chemical Bottle Wall Cabinet



#### Product Description

Our wall mounted lockable chemical storage cabinets have the capacity for 2 x 5 litre bottles of EcoWash+ wash and scale inhibitor and are available with standard key lock or combination lock.

Order Code	Description	Qty
CAB-CHEM-001	Wall cabinet 5L lock & key	1
CAB-CHEM-002	Wall cabinet 5L Combination lock	1

**Chemical Alarm Monitoring System (CAMS) takes the guess work out of when to replace the chemicals in a Panamatic Bedpan Washer, providing an audio and visual warning when the chemical bottle is empty.**

The alarm can be set to either a loud or a softer pulsed beep, and there is also the option of a night mode which emits a red flashing light as a warning so as not to disturb sleep.

### CAMS System



#### Product Description

Easy to fit the Chemical Alarm Monitoring System has the capacity to hold 2 x 5 litre bottles of EcoWash+ wash and scale inhibitor and are available with a standard key lock or combination lock.

Order Code	Description	Qty
CAB-CHEM-007	CAMS Wall cabinet 5L lock & key	1
CAB-CHEM-008	CAMS Wall cabinet 5L Combination lock	1

## Reusable Plastic Containers

All our reusable items are **lightweight strong and durable with smooth material surfaces ensuring ease of use and ease of cleaning.**

Manufactured from polypropylene or HDPE they are registered as Class 1 Medical Devices (CE Marked) and all have a shelf life of 36 months from date of manufacture providing the recommended storage conditions are adhered to.

They are designed with smooth edges and contoured shapes to ensure maximum comfort. The sturdy construction of these bedpans also helps prevent leaks or spills, providing a more sanitary experience for both patients and caregivers.

### Adult Vector Pan



#### Product Description

Reusable polypropylene adult vector bedpan suitable for both male and female use. Ideal for patients confined to bed. Lightweight, strong and durable with a secure fitting lid to prevent leaks and spills and reduces cross infection and unpleasant odours. Suitable for bariatric use, weight tested up to 250kg (39.5 stone).

Order Code	Description	UOM	Pack Qty
90-02-001	Adult Vector Pan	1	15

### Selina Super Slipper Pan



#### Product Description

Reusable polypropylene female specific slipper pan Female specific slipper urinal. Ideal for Infant use or during Maternity. Lightweight, strong and durable with a secure fitting lid to prevent leaks and spills and reduces cross infection and unpleasant odours. Can be cleaned in automatic washer disinfectors or hand washed.

Order Code	Description	UOM	Pack Qty
90-02-002	Selina Super Slipper Pan	1	25

### Favourite Fracture Pan



#### Product Description

Reusable Suitable for both male and female use, ideal for patients confined to bed. Lightweight, strong and durable with a secure fitting lid to prevent leaks and spills and reduces cross infection and unpleasant odours. Can be cleaned in automatic washer disinfectors or hand washed. Suitable for bariatric use - weight tested up to 250kg (39.5 stone).

Order Code	Description	UOM	Pack Qty
90-02-003	Favourite Fracture Pan	1	12

### Adult Male Urinal Bottle



#### Product Description

Reusable male urinal bottle. Lightweight, strong and durable with a secure fitting cap to prevent leaks and spills and reduces unpleasant odours and risk of cross infection. Graduated with clear, easy to read measurements it can be cleaned in automatic washer disinfectors or hand washed.

Order Code	Description	UOM	Pack Qty
90-02-004	Adult Male Urinal Bottle	1	45

### Standard Male Urinal Bottle



#### Product Description

Reusable standard male urinal bottle. Secure fitting cap to prevent leaks and spills and reduces unpleasant odours and risk of cross infection. Graduated with clear, easy to read measurements it can be cleaned in automatic washer disinfectors or hand washed.

Order Code	Description	UOM	Pack Qty
90-02-010	Standard Male Urinal Bottle	1	25

### Cygnet Female Urinal Bottle

**Product Description**

Reusable anatomically shaped urinal for female use. Lightweight, strong and durable with clear, easy to read graduation and easy grip handle - particularly suitable for those with limited hand movement.



Order Code	Description	UOM	Pack Qty
90-02-005	Cygnet Female Urinal Bottle	1	36

### Vomit Bowl with Handle

**Product Description**

Reusable vomit bowl with handle Ideal for patients confined to bed due to illness or surgery. Can be cleaned in automatic washer disinfectors or hand washed.



Order Code	Description	UOM	Pack Qty
90-02-008	Vomit Bowl with Handle	1	25

### Traditional Round Commode Pot

**Product Description**

Reusable traditional round commode pot with lid. Lightweight, strong and durable for use with British Standard commode chairs. Secure fitting lid to prevent leaks and spills and reduces cross infection and unpleasant odours.



Order Code	Description	UOM	Pack Qty
90-02-006	Traditional Round Commode Pot	1	25

### Wall Racks

**Product Description**

Our large wall racks hold round and square commodes, bed pans, wash bowls and bottles. This large rack will hold 5 commode & 7 urine bottle wall rack, 630mm (h) x 835mm (w) x 390mm (d).



It is wall mountable and supplied with a stainless steel drip tray.

Order Code	Description	Qty
WALLRACKASSY	Wall rack assembly 5 commode / 7 urine bottle	1

### Bedpan Holder Flat Plate

**Product Description**

Bedpan Rack designed to accommodate an individual bedpan, suitable for use with both Adult Vector Pan (90-02-001) and Selina Super Slipper Pan (90-02-002).



Order Code	Description	Qty
91-02-004	Bedpan Holder Flat Plate	1

### Wall Racks

**Product Description**

Our small wall racks hold round and square commodes, bed pans, wash bowls and bottles. This small rack holds 4 commode & 5 urine bottles, 630mm (h) x 630mm (w) x 390mm (d).



It is wall mountable and supplied with a stainless steel drip tray.

Order Code	Description	Qty
WALLRACKASSY-SMALL	Wall rack assembly 4 commode / 5 urine bottle	1



**DDC Dolphin are your  
infection control experts.**



**DDC Dolphin Ltd**  
The Fulcrum,  
Vantage Way,  
Poole,  
Dorset,  
BH12 4NU



[www.ddcdolphin.com](http://www.ddcdolphin.com)



+44 (0) 1202 731555



[info@ddcdolphin.com](mailto:info@ddcdolphin.com)



### Accreditations

